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Indications for use

Dental implant ORBITZYGOMA® is a BIONICA® dental implant system. The ORBITZYGOMA® implants are endosseous dental implants intended for surgical placement jaw in the zygomatic bone. The ORBITZYGOMA® implants are comprised of dental implant fixtures and prosthetic devices. The ORBITZYGOMA® implants provide a means for prosthetic attachment in partially and fully edentulous spans utilized delayed or immediate loading when a good primary stability is achieved in order to restore patient chewing function in the severe atrophy of the maxilla.

Contraindications

Placement of ORBITZYGOMA® dental implants may be precluded by patient's conditions that are contraindications for surgery. Absolute contraindications allegry to titanium, recent myocardial infarction, valvular prostheses, severe kidney disorders, severe liver disorders, treatment resistant diabetes, recent high-dose radiotherapy, severe chronic alcoholism, drug use, terminal illness and tumors. Relative contraindications: chemotherapy, minor kidney disorders, minor liver disorders, endome disorders, psychological and psychic disorders, lack of understanding and motivation, AIDS, HIV positive cases, prolonged use of corticosteroids, disorders of the calcium/phosphorus metabolism, erythrocytic disorders, severe nicotine poisoning. The ORBITZYGOMA® implants should not be placed in patients where the remaining zygomatic bone is too diminished or quality bone is poor to provide an adequate primary stability. The ORBITZYGOMA® implants must not be placed in patients with pathology in the maxillary sinus or in patients with anatomy defects in the maxillary sinus.

Warnings

It is recommended that clinicians, new as well experienced implant users, are updated on the BIONICA® dental implants and abutments features and their surgical applications. For detailed information about method of use of BIONICA® dental implants and abutments please consult the surgical manuals. Excessive bone loss or breakage of implant may occur when implant is loaded beyond its functional capacity. Physiological and anatomic conditions may negatively affect the performance of dental implants. This should be taken into consideration when placing dental implants in patients with the following conditions:

- poor quality bone
- poor oral hygiene
- medical conditions such as blood disorders or uncontrolled hormonal conditions
- Periodic checks for plaque on the zone of the implants through non-use recall sessions (every 4-6 months) in order to prevent the formation of periimplantitis, an undesirable event.
- The device is single use. Do not re-use and dispose after use in accordance to the regulations. Any re-use can jeopardize or seriously aggravate the health status due to risk of infections to the patient and/or the user and for the loss of functionality and efficiency of the device.
- The BIONICA® implant system has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of BIONICA® implant system in the MR environment is unknown.

Precautions

Improper technique can lead implant failure or loss of supporting bone or inadequate primary stability or other adverse effects. Through screening of implants candidates should be performed the absence of maxillar sinus pathology, an evaluation of patient capacity to maintain a correct oral hygiene, the patient motivation, the presence of parafunclional habits. Visual inspection as well as perioral radiographs is essential to determine anatomical landmarks, occlusal conditions, periodontal status and adequacy of bone. It is strongly suggested TAC dental scan to provide a real evaluation of dimensions and quality of bone.

General

The clinician must attend a specific advanced training beyond the supervision of expert implantologists about the ORBITZYGOMA® implant placement. The ORBITZYGOMA® implants are designed to reach and make use of the zygomatic bone. ORBITZYGOMA® implants can be tilted up to 45° from the occlusal plane. The implant site begins with the exposure of the maxillary lateral wall, a mucoperiosteal flap and total thickness is assessed by making a crest incision with bilateral distal vertical incision on the tuberosity areas. During the surgical exposure of the lateral wall of the jaw it is absolutely necessary to pay attention to the vital structures, including the nerves, veins and arteries. Injuries to these anatomical structures can lead to complications such as traumas of the eye, as well as extensive bleeding and dysfunctions associated with the nerves. In general the most notable risks associated with the Zyoma implants are sinusitis and fistula formations.

Radiologic examination

As with any implant patient case, a radiographic assessment is essential. As far as the ORBITZYGOMA® implant protocol is concerned, the main objectives are twofold:

1. To detect the presence of any pathology in the maxillary sinuses, bearing in mind that the thickness of the entral mucosa should not exceed 6mm.
2. To evaluate the volume of the maxillary bone surrounding the Zygomatic cavity.

Site preparation sequence

Make entrance mark with the round burr drill and continue with the drill Ø 2.8 mm until it penetrates the outer cortical layer of the zygomatic bone. Continue with the drill Ø 3.2 mm and if the bone is hard, complete the osteotomy with the drill Ø 3.6 mm.

Sterility

All implants are delivered sterile and are for single use only prior to the labeled expiration date. Do not re-sterilize. Do not use implants if the package has been damaged or previously opened. Use of non-sterile device may lead to infection of tissues or infectious diseases. For steam sterilization of the components and the Autoclave Type: Pre-Vacuum - Drying time: 20 minutes Wrap: type packaging with sealed containment device. Do not re-sterilize the same component for more than 30 times.

Potential adverse events

Possible adverse events associated with the use of dental implants may include:

1. Failure of osseointegration
2. Loss of osseointegration
3. Disconance requiring bone grafting
4. Persistent pain, numbness, paresthesia
5. Inflammation or infection
6. Excessive bone loss requiring intervention
7. Implant breakage.

Storage and Handling

Devices should be stored at room temperature. Refer to individual product labels and the surgical manual for special storage or handling conditions.

For additional information on surgical procedures, please consult the surgical manual available at www.bionica.it or request latest printed version from a BIONICA® representative.

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